



<b>ESSILOR</b>	<b>INSTRUMENTS</b> 1, rue Thomas Edison 94028 - CRETEIL Cedex	Identification <b>ILLINOIS/SYN/10</b>	Issue <b>A</b>	Page/Nbr <b>1/3</b>
----------------	---	--	-------------------	------------------------

TYPE: **SYNOPTICAL**

TITLE: **OPTEC 5000 / US Version 2**

ADDRESSEE - (Departments/Names)

**RESPONSABLE INDUSTRIEL (Stéréo Optical) A. CADOUX**  
**RESPONSABLE INDUSTRIEL (Stéréo Optical) K. AYRAULT**


	DEPT	NAME	DATE	SIGNATURE
WRITTEN BY	R&D	J. JEHAN	03/03/2004	
CHECKED BY	R&D	N. DESFARGES	03/03/2004	
APPROVED BY		NOT APPLICABLE		

NAME and POST of MANAGER  
Laurence. LE REOUR

Issue	Date	Subject of modification
A	03/03/2004	CREATION OF DOCUMENT.

SHEET (Page)	HISTORY	Date	Rev.
	CREATION OF DOCUMENT	3/03/2004	A

<b>ESSILOR INSTRUMENTS</b> 64Bis, Avenue Aubert 94300 VINCENNES TEL: 01.72.94.71.00 FAX: 01.72.94.70.24 ILLINOIS			
<b>SYNOPTICAL OPTEC 5000 / US - V2</b>			
Size B	Document Number	SYNOPTICAL OPTEC 5000US V2A.DSN	Rev A
Date: Wednesday, March 03, 2004	Sheet	1 of 3	EIC/R&D/FDQ/07-C



32110 CORD  
POWER SUPPLY  
HOSPITAL  
GRADE

**POWER SUPPLY**  
24VDC/1A  
32083 POWER SUPPLY

32101 CABLE  
SUBASSEMBLY  
POWER JACK  
IL1R24

OPTION NC  
JP1 P21  
**52046 PCB  
Control  
Panel assm.**  
(IL2R01) P1  
OPTION NC  
OPTION NC

IL1R11  
32089 CABLE LIGHT READY INDEXATION

32099 CABLE  
POWER SWITCH  
ASSEMBLY  
IL1R21


**52053 PCB  
Indexation  
assm.**  
Top Test  
(IL1R05)

**ESSILOR INSTRUMENTS**  
64Bis, Avenue Aubert 94300 VINCENNES  
TEL: 01.72.94.71.00 FAX: 01.72.94.70.24  
ILLINOIS

Title <b>SYNOPTICAL OPTEC 5000 US / - V2</b>		
Size B	Document Number SYNOPTICAL OPTEC 5000US V2.A.DSN	Rev A
Date: Wednesday, March 03, 2004	Sheet 2 of 3	EIC/R&D/FDQ/07-C



QUICK-CONNECT 6.35mm x 0.8  
 ILLR20  
 32098 CABLE READY SWITCH



32110 CORD  
 POWER SUPPLY  
 HOSPITAL  
 GRADE

**POWER SUPPLY FOR  
 MEDICAL USE**  
 AC INPUT DC OUPUT  
 100-240V 24V  
 0.6-0.3A  
 50-60HZ

32083 POWER SUPPLY

32101 CABLE  
 SUBASSEMBLY  
 POWER JACK  
 ILLR24

**52046 PCB  
 Control  
 Panel assm.**  
 (IL2R01)

**52049 PCB Target  
 Light assm.**  
 (IL1R03)

FAR RIGHT

**52049 PCB Target  
 Light assm.**  
 (IL1R03)

FAR LEFT

**52049 PCB Target  
 Light assm.**  
 (IL1R03)

NEAR RIGHT

**52049 PCB Target  
 Light assm.**  
 (IL1R03)

NEAR LEFT

**52053 PCB  
 Indexation assm.  
 Top Test**  
 (IL1R05)

**ESSILOR INSTRUMENTS**  
 64Bis, Avenue Aubert 94300 VINCENNES  
 TEL: 01.72.94.71.00 FAX: 01.72.94.70.24  
 ILLINOIS

Title <b>SYNOPTICAL OPTEC 5000 / US - V2détail</b>		
Size B	Document Number SYNOPTICAL OPTEC 5000US V2A.DSN	Rev A
Date: Wednesday, March 03, 2004	Sheet 3 of 3	EIC/R&D/FDQ07-C

POWER SWITCH  
 ON/OFF  
 CONTROL  
 POWER  
 LED  
 GREEN

32099 CABLE POWER  
 SWITCH ASSEMBLY  
 ILLR21