

RECORD FORM - F.A.C.T.® - FUNCTIONAL VISION ANALYZER™

NAME: _____
 ADDRESS: _____

 ADMINISTRATOR: _____

DATE: _____
 AGE: _____
 WEARING GLASSES: Y N
 CONTACT LENSES: Y N

DIAL AT 01: FAR ACUITY MONOCULAR

LINE/ ACUITY	LEFT	RIGHT
5 20/160	DRCHV	RVCSH
6 20/125	CKNRD	CKDZR
7 20/100	SHZDO	OVRHK
8 20/80	RODVC	NRVKO
9 20/63	KRHSD	KSNDK
10 20/50	COSZH	VHCRD
11 20/40	ZCVOR	DSRKH
12 20/32	CRDVH	KRSND
13 20/25	DCVHS	SZVHO
14 20/20	KVSCR	HRCNS
15 20/16	OCNKD	ZCVNO
16 20/12.5	DKCVZ	OKZHC

DIAL AT 02: FAR ACUITY BINOCULAR

LINE/ ACUITY	BOTH EYES
5 20/160	SCNZV
6 20/125	CSHDN
7 20/100	ONKCH
8 20/80	CVZHO
9 20/63	VCHON
10 20/50	RDCZK
11 20/40	HOSDR
12 20/32	RSOVH
13 20/25	NOKDR
14 20/20	ZHSOK
15 20/16	CDKVH
16 20/12.5	HKDCO

DIAL AT 03: NEAR ACUITY MONOCULAR

LINE/ ACUITY	LEFT EYE	RIGHT EYE
5 20/160	RDKVC	SCNZV
6 20/125	VOSDR	CSHDN
7 20/100	OZCRS	ONKCH
8 20/80	COKRN	CVZHO
9 20/63	ZHSKO	VCHON
10 20/50	CVDHK	RDCZK
11 20/40	RHODS	HOSDR
12 20/32	ONVCZ	RSOVH
13 20/25	HDRVC	NOKDR
14 20/20	KRNSD	ZHSOK
15 20/16	CZDVK	CDKVH
16 20/12.5	HDKCN	HKDCO

DIAL AT 04: NEAR ACUITY BINOCULAR

LINE/ ACUITY	BOTH EYES
5 20/160	RDKVC
6 20/125	VOSDR
7 20/100	OZCRS
8 20/80	COKRN
9 20/63	ZHSKO
10 20/50	CVDHK
11 20/40	RHODS
12 20/32	ONVCZ
13 20/25	HDRVC
14 20/20	KRNSD
15 20/16	CZDVK
16 20/12.5	HDKCN

DIAL AT 05 through 09:

Functional Acuity Contrast Test
 Charts on Reverse Side

DIAL AT 07: FIVE F.A.C.T. SLIDE - LINE C

CATARACT FUNCTIONAL DISABILITY TEST

	Pass	Fail
Left	<input type="checkbox"/>	<input type="checkbox"/>
Right	<input type="checkbox"/>	<input type="checkbox"/>

NIGHT WITHOUT GLARE Target = 3cd/m ²	NIGHT WITH GLARE Target = 3cd/m ² Glare Test 1 = 1 Lux	DAY WITHOUT GLARE Target = 85cd/m ²	DAY WITH GLARE Target = 85cd/m ² Glare Test 1 = 10 Lux
9 ●	9 ●	9 ●	9 ●
8 ●	8 ●	8 ●	8 ●
7 ●	7 ●	7 ●	7 ●
6 ●	6 ●	6 ●	6 ●
5 ●	5 ●	5 ●	5 ●
4 ●	4 ●	4 ●	4 ●
3 ●	3 ●	3 ●	3 ●
2 ●	2 ●	2 ●	2 ●
1 ●	1 ●	1 ●	1 ●
C	C	C	C

■ If the subject falls in any gray shaded area, a full contrast sensitivity test should be given.

DIAL AT 10: DEPTH PERCEPTION - FAR

1-B	2-L	3-B	4-T	5-T	6-L	7-R	8-L	9-R
400	200	100	70	50	40	30	25	20

DIAL AT 11: COLOR PERCEPTION - FAR

A-12	B-5	C-2.6	D-6	E-16	F-BLANK
PASS	_____				FAIL

DIAL AT 12: FAR ACUITY MONOCULAR

LINE ACUITY	LEFT EYE	RIGHT EYE
5 20/160	DRCHV	RVCSH
6 20/125	CKNRD	CKDZR
7 20/100	SHZDO	OVRHK
8 20/80	RODVC	NRVKO
9 20/63	KRHSD	KSNDK
10 20/50	COSZH	VHCRD
11 20/40	ZCVOR	DSRKH
12 20/32	CRDVH	KRSND
13 20/25	DCVHS	SZVHO
14 20/20	KVSCR	HRCNS
15 20/16	OCNKD	ZCVNO
16 20/12.5	DKCVZ	OKZHC




STEREO OPTICAL

3539 N. Kerton Avenue
 Chicago, IL 60641-38798 USA
 Phone: 1.773.777.2869 (in Illinois) or 1.800.344.9500 (USA & Canada)
 Fax: 1.773.777.4985
 Email: sales@stereooptical.com
www.stereooptical.com or www.functionalvisionanalyzer.com


© Copyright Stereo Optical Co., Inc 2006;
 Functional Vision Analyzer is a trademark of Stereo Optical Co., Inc.
 F.A.C.T. is a registered trademark of Stereo Optical Co., Inc.

NIGHT TESTING WITHOUT GLARE




NIGHT TESTING WITH GLARE

- Glare Test 1 (1 Lux)
- Glare Test 2 (28 Lux)



DAY TESTING WITHOUT GLARE



DAY TESTING WITH GLARE

- Glare Test 1 (10 Lux)
- Glare Test 2 (135 Lux)

