


Stereo Optical Slide Library




2000-001
Far Acuity Binocular Landolt Rings
(20/200,20/100,20/70,20/50,20/40,20/35,20/30,
20/25,20/22,20/20,20/18,20/17,20/15,20/13)



2000-002
Near Acuity Binocular Landolt Rings
(20/200,20/100,20/70,20/50,20/40,20/35,20/30,
20/25,20/22,20/20,20/18,20/17,20/15,20/13)



2000-003
Far Acuity Binocular Letter
(20/200,20/100,20/70,20/50,20/40,20/30,20/20)



2000-004
Far Acuity Tumbling E Monocular
(20/200,20/100,20/70,20/50,20/40,20/30,20/20)



2000-005
Far Acuity Monocular Left Eye Landolt Rings
(20/200,20/100,20/70,20/50,20/40,20/35,20/30,
20/25,20/22,20/20,20/18,20/17,20/15,20/13)



2000-006
Near Acuity Monocular Left Eye Landolt Rings
(20/200,20/100,20/70,20/50,20/40,20/35,20/30,
20/25,20/22,20/20,20/18,20/17,20/15,20/13)


Slide backgrounds altered for clarity of targets




2000-007
Near Acuity Binocular Letter
(20/100,20/70,20/50,20/40,20/30,20/20)



2000-008
Far Acuity Monocular Right Eye Landolt Rings
(20/200,20/100,20/70,20/50,20/40,20/35,20/30,
20/25,20/22,20/20,20/18,20/17,20/15,20/13)



2000-009
Near Acuity Binocular Landolt Rings
(20/200,20/100,20/70,20/50,20/40,20/35,20/30,
20/25,20/22,20/20,20/18,20/17,20/15,20/13)



2000-010
Far Color Perception Binocular
(Pseudo Ishihara)




2000-011
Far Demonstration Mountains Binocular
(Fixation for peripheral test)



2000-012
Far Lateral Phoria Binocular
(1 Diopter Increments)

Slide backgrounds altered for clarity of targets



2000-013
Near Lateral Phoria Binocular
(1-1/2 Diopter Increments)



2000-014
Far Acuity Michigan Preschool Binocular
20/70



2000-015
Far Acuity Michigan Preschool Monocular
20/30



2000-016
Far Acuity Michigan Preschool Monocular
20/30




2000-017
Far Acuity Michigan Preschool Monocular
20/30



2000-018
Far Acuity Michigan Preschool Monocular
20/30


Slide backgrounds altered for clarity of targets



2000-019
Far Pediatric Color Perception Binocular
Tumbling "E"




2000-020
Far Muscle Balance
(Combination Lateral & Vertical Phoria)



2000-021
Far Multiple Phoria Binocular
(2.3 Diopter Horizontal
3.6 Diopter Vertical)



2000-022
Near Muscle Balance
(Combination Lateral & Vertical Phoria)



2000-023
Near Multiple Phoria Binocular
(2.7 Diopter Horizontal 3.5 Diopter Vertical)



2000-024
Far Stereo Depth Perception
(400,200,100,70,50,40,30,25,20 Seconds of Arc)

Slide backgrounds altered for clarity of targets



2000-025
Far Vertical Phoria Binocular
(1/2 Diopter Increments)



2000-026
Far Acuity HOTV Binocular
20/70



2000-027
Far Acuity HOTV Monocular
20/40



2000-028
Far Acuity HOTV Monocular
20/40



2000-029
Far Acuity HOTV Monocular
20/40




2000-030
Far Acuity HOTV Monocular
20/30

Slide backgrounds altered for clarity of targets



3000-038
50% Far Acuity Monocular
C.A.T.® Contrast Sensitivity
 (20/160,20/125,20/100,20/80,20/63,20/50,20/40,
 20/32,20/25,20/20,20/16, 20/12.5)
(Slides 037 to 041 sold as set only)



3000-039
25% Far Acuity Monocular
C.A.T.® Contrast Sensitivity
 (20/160,20/125,20/100,20/80,20/63,20/50,20/40,
 20/32,20/25,20/20,20/16, 20/12.5)
(Slides 037 to 041 sold as set only)



3000-040
12.5% Far Acuity Monocular
C.A.T.® Contrast Sensitivity
 (20/160,20/125,20/100,20/80,20/63,20/50,20/40,
 20/32,20/25,20/20,20/16, 20/12.5)
(Slides 037 to 041 sold as set only)



3000-041
6.25% Far Acuity Monocular
C.A.T.® Contrast Sensitivity
 (20/160,20/125,20/100,20/80,20/63,20/50,20/40,
 20/32,20/25,20/20,20/16, 20/12.5)
(Slides 037 to 041 sold as set only)



2000-042
Far Acuity Monocular
 (20/160,20/125,20/100,20/80,20/63,20/50,20/40,
 20/32,20/25,20/20,20/16, 20/12.5)




2000-043
Far Astigmatic Clock Binocular


Slide backgrounds altered for clarity of targets




2000-050
Far Acuity Binocular Letters
(20/100,20/75,20/60,20/50,20/40,20/20)




2000-051
Near Acuity Binocular Letters
(20/200,20/100,20/70,20/50,20/40,20/30,20/20)




2000-052
Far Acuity Binocular Letters
(20/100,20/70,20/50,20/40,20/30,20/25,20/20,20/15)



2000-054
Far Acuity Binocular Letters
(20/100,20/70,20/50,20/40,20/30,20/20)




2000-055
Far Acuity Binocular Letter
20/40 all lines




2000-056
Far Acuity Binocular Letters
Lines 1- 4 20/40
Lines 5-6 20/100


Slide backgrounds altered for clarity of targets



2000-065
Far Acuity Binocular Numbers
(20/200,20/100,20/70,20/50,20/40,20/30,20/20)



2000-066
Near Acuity Binocular Letters
Top Row (20/200,20/100)
(20/70,20/60,20/50,20/40,20/30,20/25,20/20)



2000-067
Far Acuity Binocular Numbers
Top Row (20/200,20/100)
(20/70,20/60,20/50,20/40,20/30,20/25,20/20)



2000-068
Far Acuity Monocular Numbers
Top Row (20/200,20/100,20/70)
(20/50,20/40,20/30,20/25,20/20,20/15)




2000-070
Far Acuity Monocular Tumbling "E"
Top Row (20/200,20/100)
(20/70,20/50,20/40,20/30,20/20)




2000-071
Far Acuity Binocular Tumbling "E"
(20/50,20/40,20/30,20/20)

Slide backgrounds altered for clarity of targets



2000-072
 Near Acuity Allen Test
 Binocular 20/100
 Monocular (20/50,20/40,20/30)



2000-073
 Near Acuity Binocular Letters
 (20/100,20/70,20/50,20/40,20/30,20/20)



2000-074
 Far Acuity Monocular Tumbling "E"
 Top Row 20/200,20/100,20/70,20/50,20/40,20/30,20/20
 (20/50,20/40,20/30,20/20)



2000-075
 Near Near Point Reading
 (.37, .50, .62, .75 M)



2000-077
 Far Point Color Perception (Pseudo Ishihara)
 Far Binocular Acuity Letter 20/40
 Far Binocular Acuity Landolt "C" 20/40
 Far Stereo Depth Perception (24 & 86 Seconds of Arc)



2000-077C
 Far Point color Perception (Pseudo Ishihara)
 Far Stereo Depth Perception (24 & 86 Seconds of Arc)

Slide backgrounds altered for clarity of targets




2000-079
 Far Acuity Binocular Letters 20/40
 Road Sign Recognition
 Color Recognition




2000-080
 Far Acuity Binocular Letters 20/40
 Far Road Sign Recognition
 Far Color Perception (Pseudo Ishihara)



2000-081
 Near Acuity Binocular Letters 20/40
 Far Road Sign Recognition
 Far Color Perception (Pseudo Ishihara)



2000-082
 Far Color Perception Binocular
 (Pseudo Ishihara)



2000-084
 Far Color Recognition Binocular
 (Red, Green, Amber)




2000-087
 Far Fusion Binocular

Slide backgrounds altered for clarity of targets



2000-088
Near Fusion Binocular



2000-089
Far Acuity Binocular Landolt "C"
(20/100,20/75,20/60,20/50,20/40,20/20)



2000-090
Far Acuity Binocular Michigan Preschool
20/75



2000-091
Near Acuity Binocular Michigan Preschool
20/70



2000-092
Far Acuity Monocular Michigan Preschool
20/20



2000-093
Far Acuity Monocular Michigan Preschool
20/20


Slide backgrounds altered for clarity of targets




2000-094
Far Acuity Binocular Michigan Preschool
20/20



2000-095
Far Acuity Monocular Michigan Preschool
20/20



2000-096
Far Acuity Monocular Michigan Preschool
20/40



2000-097
Far Acuity Monocular Michigan Preschool
20/40



2000-098
Far Acuity Monocular Michigan Preschool
20/40



2000-099
Far Acuity Monocular Michigan Preschool
20/40

Slide backgrounds altered for clarity of targets



2000-100
Far Acuity Monocular Michigan Preschool
20/50



2000-101
Far Acuity Monocular Michigan Preschool
20/50



2000-102
Far Acuity Monocular Michigan Preschool
20/50



2000-103
Far Acuity Monocular Michigan Preschool
20/50



2000-104
Far Acuity Monocular Michigan Preschool
20/70



2000-105
Far Acuity Monocular Michigan Preschool
20/70

Slide backgrounds altered for clarity of targets



2000-106
Far Acuity Binocular Michigan Preschool
20/70



2000-107
Far Acuity Monocular Michigan Preschool
20/70



2000-108
Far Acuity Monocular Michigan Preschool
20/100



2000-109
Far Acuity Monocular Michigan Preschool
20/100



2000-110
Far Acuity Binocular Michigan Preschool
20/100



2000-111
Far Acuity Monocular Michigan Preschool
20/100

Slide backgrounds altered for clarity of targets



2000-115
Far Phoria Binocular
(Picnic Table with red apple)




2000-116
Far Phoria Binocular




2000-117
Near Phoria Binocular
(Picnic Table with red apple)



2000-118
Near Phoria Binocular



2000-120
Far Stereo Depth Perception
(400,200,100,70,50,40,30,25,20 Seconds of Arc)




2000-125
Near Acuity Binocular Letters
Top Row (20/200,20/100,20/70)
(20/50,20/40,20/30,20/25,20/20,20/15)

Slide backgrounds altered for clarity of targets



2000-126
Far Acuity Monocular Letters
Top Row (20/200,20/100,20/70)
(20/50,20/40,20/30,20/25,20/20,20/15)



2000-128
Near Acuity Binocular Letters
Top Row (20/200,20/100,20/70)
(20/50,20/40,20/30,20/25,20/20,20/15)



2000-129
Far Acuity Monocular Letters
Top Row (20/200,20/100,20/70)
(20/50,20/40,20/30,20/25,20/20,20/15)



2000-132
Far Road Sign Recognition
Depth Perception Binocular
(620, 345 & 155 Seconds of Arc)




2000-137
Far Acuity Low Vision "C"
(20/700,20/600,20/500,20/400,20/300,20/200,20/100)

Slide backgrounds altered for clarity of targets



2000-138
Far Acuity Low Vision Tumbling “E”
(20/800,20/600,20/400)



2000-145
Near Vertical Phoria Binocular
(1/2 Prism Diopter)




2000-147
Far Demonstration Highway Binocular
(Fixation for peripheral test)



2000-148
Far Road Sign Recognition &
Stereo Depth Perception Binocular
(620,323 & 145 Seconds of Arc)




2000-149
Far Color Recognition
(Red, Yellow, Green, Blue)




3000-171
Far F.A.C.T.[®] Contrast Sensitivity Monocular
1.5 Cycles Per Degree
(Slides 171 to 175 sold as set only)


Slide backgrounds altered for clarity of targets




3000-172
Far F.A.C.T.[®] Contrast Sensitivity Monocular
3 Cycles Per Degree
(Slides 171 to 175 sold as set only)




3000-173
Far F.A.C.T.[®] Contrast Sensitivity Monocular
6 Cycles Per Degree
(Slides 171 to 175 sold as set only)




3000-174
Far F.A.C.T.[®] Contrast Sensitivity Monocular
12 Cycles Per Degree
(Slides 171 to 175 sold as set only)



3000-175
Far F.A.C.T.[®] Contrast Sensitivity Monocular
18 Cycles Per Degree
(Slides 171 to 175 sold as set only)




3000-176
Near F.A.C.T.[®] Contrast Sensitivity Monocular
1.5 Cycles Per Degree
(Slides 176 to 178 sold as set only)




3000-177
Near F.A.C.T.[®] Contrast Sensitivity Monocular
3 & 6 Cycles Per Degree
(Slides 176 to 178 sold as set only)

Slide backgrounds altered for clarity of targets



3000-178
 Near F.A.C.T.[®] Contrast Sensitivity Monocular
 12 & 18 Cycles Per Degree
(Slides 176 to 178 sold as set only)



2000-182 Far Visual Acuity Binocular HOTV
 (20/100,20/40,20/32,20/25,20/20)



2000-183 Far Visual Acuity Monocular HOTV
 (20/100,20/40,20/32,20/25,20/20)



2000-184 Near Visual Acuity Numbers
 (20/40)



2000-187 Near Visual Acuity Binocular HOTV
 (20/100,20/40,20/32,20/25,20/20)



2000-188 Near Visual Acuity Monocular HOTV
 (20/100,20/40,20/32,20/25,20/20)

Slide backgrounds altered for clarity of targets



2000-253
 Far Vertical Phoria Binocular
 (1/2 Diopter Increments)



2000-254
 Far Lateral Phoria Binocular
 (1/2 Diopter Increments)



2000-255
 Far Acuity Binocular Checkerboard
 (20/200,20/100,20/67,20/50,20/40,20/33,20/29,
 20/25,20/22,20/20,20/18,20/17)



2000-256
 Far Acuity Monocular Right Eye Checkerboard
 (20/200,20/100,20/67,20/50,20/40,20/33,20/29,
 20/25,20/22,20/20,20/18,20/17)



2000-257
 Far Acuity Monocular Left Eye Checkerboard
 (20/200,20/100,20/67,20/50,20/40,20/33,20/29,
 20/25,20/22,20/20,20/18,20/17)



2000-258
 Far Stereo Depth Perception
 (362,83,43,32,27,19,13,12,9.7 Seconds of Arc)

Slide backgrounds altered for clarity of targets




2000-259
Far Color Perception
(Pseudo Ishihara)



2000-260
Near Acuity Binocular Checkerboard
(20/200,20/100,20/67,20/50,20/40,20/33,20/29,
20/25,20/22,20/20,20/18,20/17)



2000-261
Near Acuity Monocular Right Eye Checkerboard
(20/200,20/100,20/67,20/50,20/40,20/33,20/29,
20/25,20/22,20/20,20/18,20/17)



2000-262
Far Acuity Monocular Left Eye Checkerboard
(20/200,20/100,20/67,20/50,20/40,20/33,20/29,
20/25,20/22,20/20,20/18,20/17)



2000-263
Near Vertical Phoria Binocular
(1/2 Diopter Increments)



2000-264
Near Lateral Phoria Binocular
(1 Diopter Increments)

Slide backgrounds altered for clarity of targets




2000-265
Far Acuity Binocular Tumbling "E"
(20/200,20/100,20/50,20/40,20/30,20/20)



2000-266
Far Phoria Binocular




2000-267
Near Phoria Binocular



2000-268
Far Acuity Binocular Letter
(20/200,20/100,20/70,20/50,20/40,20/30,20/20)


Slide backgrounds altered for clarity of targets




2000-185
 NEAR Acuity Monocular Letters
 (20/160, 20/125, 20/100, 20/80, 20/63, 20/50, 20/40, 20/32, 20/25, 20/20, 20/16, 20/12.5)




2000-189
 NEAR Acuity Binocular Letters
 (20/160, 20/125, 20/100, 20/80, 20/63, 20/50, 20/40, 20/32, 20/25, 20/20, 20/16, 20/12.5)




2000-191
 Far Binocular Contrast Letters
 25% ETDRS
 (20/160, 20/125, 20/100, 20/80, 20/63, 20/50, 20/40, 20/32, 20/25, 20/20, 20/16, 20/12.5)



2000-192
 Far Binocular Contrast Letters
 12.5% ETDRS
 (20/160, 20/125, 20/100, 20/80, 20/63, 20/50, 20/40, 20/32, 20/25, 20/20, 20/16, 20/12.5)



2000-193
 Far Binocular Contrast Letters
 50% ETDRS
 (20/160, 20/125, 20/100, 20/80, 20/63, 20/50, 20/40, 20/32, 20/25, 20/20, 20/16, 20/12.5)



2000-195
 NEAR Acuity Binocular Letters
 100% ETDRS
 (20/160, 20/125, 20/100, 20/80, 20/64, 20/50, 20/40, 20/32, 20/25, 20/20, 20/16, 20/12.5)

Slide background altered for clarity of targets