

SPECIAL STEREO OPTICAL RECORD FORM FOR THE OPTEC VISION TESTER

NAME: _____
 ADDRESS: _____
 ADDRESS: _____
 ADMINISTRATOR: _____

DATE: _____
 AGE: _____
 WEARING GLASSES: Y N
 CONTACT LENSES: Y N

DIAL AT 01:	LINE/ACUITY	LEFT	RIGHT
<i>FAR</i>	5 20/160	DRCHV	RVCSH
<i>ACUITY</i>	6 20/125	CKNRD	CKDZR
<i>MONOCULAR</i>	7 20/100	SHZDO	OVRHK
	8 20/80	RODVC	NRVKO
	9 20/63	KRHSD	KSNDK
	10 20/50	COSZH	VHCRD
	11 20/40	ZCVOR	DSRKH
	12 20/32	CRDVH	KRSND
	13 20/25	DCVHS	SZVHO
	14 20/20	KVSCR	HRCNS
	15 20/16	OCNKD	ZCVNO
	16 20/12.5	DKCVZ	OKZHC

DIAL AT 02:	LINE/ACUITY	BOTH
<i>FAR</i>	5 20/160	SCNZV
<i>ACUITY</i>	6 20/125	CSHDN
<i>BINOCULAR</i>	7 20/100	ONKCH
	8 20/80	CVZHO
	9 20/63	VCHON
	10 20/50	RDCZK
	11 20/40	HOSDR
	12 20/32	RSOVH
	13 20/25	NOKDR
	14 20/20	ZHSOK
	15 20/16	CDKVH
	16 20/12.5	HKDCO

DIAL AT 03:	LINE/ACUITY	LEFT	RIGHT
<i>NEAR</i>	5 20/160	RDKVC	SCNZV
<i>ACUITY</i>	6 20/125	VOSDR	CSHDN
<i>MONOCULAR</i>	7 20/100	OZCRS	ONKCH
	8 20/80	COKRN	CVZHO
	9 20/63	ZHSKO	VCHON
	10 20/50	CVDHK	RDCZK
	11 20/40	RHODS	HOSDR
	12 20/32	ONVCZ	RSOVH
	13 20/25	HDRVC	NOKDR
	14 20/20	KRNSD	ZHSOK
	15 20/16	CZDVK	CDKVH
	16 20/12.5	HDKCN	HKDCO

DIAL AT 04:	LINE/ACUITY	BOTH
<i>NEAR</i>	5 20/160	RDKVC
<i>ACUITY</i>	6 20/125	VOSDR
<i>BINOCULAR</i>	7 20/100	OZCRS
	8 20/80	COKRN
	9 20/63	ZHSKO
	10 20/50	CVDHK
	11 20/40	RHODS
	12 20/32	ONVCZ
	13 20/25	HDRVC
	14 20/20	KRNSD
	15 20/16	CZDVK
	16 20/12.5	HDKCN

DIAL AT 05 thru 09: **Functional Acuity Contrast Test**
 Charts on Reverse Side

DIAL AT 10: **DEPTH PERCEPTION - FAR** 1-B 2-L 3-B 4-T 5-T 6-L 7-R 8-L 9-R
 SECONDS OF ARC 400 200 100 70 50 40 30 25 20

DIAL AT 11: **COLOR PERCEPTION - FAR** A-12 B-5 C-26 D-6 E-16 F-BLANK


DIAL AT 12:	LINE/ACUITY	LEFT	RIGHT
<i>FAR</i>	5 20/160	DRCHV	RVCSH
<i>ACUITY</i>	6 20/125	CKNRD	CKDZR
<i>MONOCULAR</i>	7 20/100	SHZDO	OVRHK
	8 20/80	RODVC	NRVKO
	9 20/63	KRHSD	KSNDK
	10 20/50	COSZH	VHCRD
	11 20/40	ZCVOR	DSRKH
	12 20/32	CRDVH	KRSND
	13 20/25	DCVHS	SZVHO
	14 20/20	KVSCR	HRCNS
	15 20/16	OCNKD	ZCVNO
	16 20/12.5	DKCVZ	OKZHC




STEREO OPTICAL Co., Inc.

3539 N Kenton Avenue
 Chicago, IL 60641-3879
 800-344-9500
 In IL 773-777-2869
 Fax 773-777-4985
 email:sales@stereooptical.com
 www.stereooptical.com


NIGHT TESTING WITHOUT GLARE



NIGHT TESTING WITH GLARE



DAY TESTING WITHOUT GLARE



DAY TESTING WITH GLARE

