




International: For complete occupational screening using accurate and easily understood targets to help determine effective vision. Includes Job Standards Booklet to easily identify individuals requiring further testing for all occupations. Intermediate Lenses available for VDT operators, which is the fastest growing population in the industry today. Now test Distance, Intermediate, and Near Point Vision.




Slide #1 (2000-255) Far Point
Acuity Both Eyes Binocular
(20/200-20/17)




Slide #7(2000-254) Far Point
Lateral Phoria (1 Diopter Increments)




Slide #2 (2000-256) Far Point
Acuity Right Eye Monocular
(20/200-20/17)



Slide #8 (2000-260) Near Point
Acuity Both Eyes Binocular
(20/200-20/17)




Slide #3 (2000-257) Far Point
Acuity Left Eye Monocular
(20/200-20/17)




Slide #9 (2000-261) Near Point
Acuity Right Eye Monocular
(20/200-20/17)




Slide #4 (2000-258) Far Point
Stereo Depth Perception
(362-9.7 Seconds of Arc)



Slide #10 (2000-262) Near Point
Acuity Left Eye Monocular
(20/200-20/17)



Slide #5 (2000-259) Far Point
Color Perception (Pseudo Ishihara)



Slide #11 (2000-264) Near Point
Lateral Phoria (1 Diopter Increments)



Slide #6 (2000-253) Far Point
Vertical Phoria (1/2 Diopter Increments)



Slide #12 (2000-263) Near Point
Vertical Phoria (1/2 Diopter Increments)