

# STEREO OPTICAL PRESSCREENING VISION TESTER RECORD FORM

NAME: \_\_\_\_\_  
 ADDRESS: \_\_\_\_\_  
 ADDRESS: \_\_\_\_\_  
 ADMINISTRATOR: \_\_\_\_\_

DATE: \_\_\_\_\_  
 AGE: \_\_\_\_\_  
 WEARING GLASSES: Y N  
 CONTACT LENSES: Y N

**PERIPHERAL TEST** L = 85 70 55 45 NASAL 45 R = 85 70 55

**DIAL AT 01: VISUAL ACUITY - FAR** Score: 20/\_\_\_

LINE/ACUITY	LEFT EYE	BOTH EYES	RIGHT EYE
1 20/200	ZN	RO	HK
2 20/100	RKS	HNC	ZOD
3 20/70	HCDV	SKZO	RNDS
4 20/50	ZROD	NSCH	VZKN
5 20/40	KHSC	OZNR	DNVC
6 20/30	ONRZV	DKHCS	KDSO
7 20/20	SDCHN	VRZKO	HSNRD

**DIAL AT 03: VISUAL ACUITY TESTS - NEAR** Score: 20/\_\_\_

LINE/ACUITY	LEFT EYE	BOTH EYES	RIGHT EYE
1 20/100	SVC	NRK	HZO
2 20/70	RNZH	DOKV	CSZN
3 20/50	CKVD	SNZR	DOHC
4 20/40	VHRN	ODSK	NZCS
5 20/30	HSKRC	NZDOV	ZSHNK
6 20/20	ZONVR	HCSKD	VKCDS

**DIAL AT 02: COLOR PERCEPTION - FAR** A-1 2 B-5 C-2 6 D-6 E-1 6 F-BLANK PASS \_\_\_ FAIL \_\_\_

**DIAL AT 04: LATERAL PHORIA - FAR** 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15  
ESOPHORIA ORTHO EXOPHORIA

**DIAL AT 05: VERTICAL PHORIA - FAR** 1 2 3 4 5 6 7  
LEFT HYPERPHORIA ORTHO RIGHT HYPERPHORIA

**DIAL AT 06: LATERAL PHORIA - NEAR** 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15  
ESOPHORIA ORTHO EXOPHORIA

**DIAL AT 07: DEPTH PERCEPTION - FAR** 1-B 2-L 3-B 4-T 5-T 6-L 7-R 8-L 9-R  
SECONDS OF ARC 400 200 100 70 50 40 30 25 20

**DIAL AT 08: FUSION - FAR** 3 cubes BOTH eyes = PASS Other than 3 cubes = FAIL

<b>DIAL AT 09: LEFT EYE score: 20/___</b>		Right Occluder OFF		<b>RIGHT EYE score: 20/___</b>		Left Occluder OFF									
<b>TUMBLING "E" FAR</b>	<b>Pilot</b>	<u>20</u>	<u>20</u>	<u>20</u>	<u>20</u>	<u>20</u>	<u>20</u>								
<b>LINE/ACUITY</b>	<b>COLOR</b>	200	100	70	50	40	30								
20/50	BLUE	L	T	B	L	R	L	T	R	B	R	L	L	T	B
20/40	RED	R	T	R	L	L	R	L	R	T	L	B	T		
20/30	AMBER	B	R	B	L	L	R	T	T	L	B	T	R		
20/20	GREEN	L	T	T	L	R	B	L	B	R	R	L	T		
		R	L	T	L	B	T	R	T	R	L	T	B		

**DIAL AT 10: TUMBLING "E" COLOR PERCEPTION - FAR** 1-R 2-L 3-B 4-T 5-B 6-L 7-T 8-R  
R=Right L=Left T=Top B=Bottom

**DIAL AT 11: MUSCLE BALANCE - FAR** Ball IN Box = PASS Ball OUT Box = FAIL

**DIAL AT 12: ALLEN PRESCHOOL - FAR**  
Both Eyes 20/100 : Car, Phone, Cake, Horse

LINE/ACUITY	LEFT EYE	RIGHT EYE
A 20/50 (BLUE)	Horse, Cake, Phone, Car	Cake, Horse, Car, Phone
B 20/40 (RED)	Cake, Car, Horse, Phone	Phone, Cake, Horse, Car
C 20/30 (YELLOW)	Phone, Horse, Car, Cake	Horse, Car, Phone, Cake



**STEREO OPTICAL CO., INC.**  
 8623 W. Bryn Mawr Avenue, Suite 502  
 Chicago, IL 60631-3501 U.S.A.  
 800-344-9500 or 773-867-0380  
 Fax 773-867-0388  
 email:sales@stereooptical.com  
 www.stereooptical.com