

Optec[®] 900 Color Vision Tester

NSN: 6540-01-452-8207

Meets the Color Perception Requirements stated in Articles 15-40 (1) (K), 15-70 (2) (G), and 15-69 (2) of "Manual of the Medical Department" (MANMED) issued by the U.S. Navy.



*10% Smaller and Weighs 35% Less
Than Previous ColorVision Tester*

- High Impact Fiber Filled Water Resistant Housing
- Self Contained Electrical Cord for Storage
- Plastic Laminated Scoring Template and Test Instructions
- Carrying Case Included
- Built-in Electronic Timer with Over-ride Switch
(Micro-processor Controlled)




**STEREO
OPTICAL**

Optec is a registered trademark of Stereo Optical Company, Inc.

www.stereooptical.com / sales@stereooptical.com / 1.773.867.0380 / 1.800.344.9500 / Fax: 1.773.867.0388

Optec[®] 900 Color Vision Tester

NSN: 6540-01-452-8207



Optec [®] 900 ColorVision Tester Includes:	Specifications:
<ul style="list-style-type: none"> • 1 ColorVision Tester with Carrying Case • Laminated Instructions • Scoring Template • Dust Cover • Service Manual 	<ul style="list-style-type: none"> • Delivery 30 Days After Receipt of Order • Underwriters Laboratory 544 and CUL Listing • 120 VAC-50/60 Hz-40 WATTS • One Year Warranty



**STEREO
OPTICAL**

Optec is a registered trademark of Stereo Optical Company, Inc.

www.stereooptical.com / sales@stereooptical.com / 1.773.867.0380 / 1.800.344.9500 / Fax: 1.773.867.0388