

STEREO OPTICAL MEDICAL RECORD FORM FOR THE OPTEC® VISION TESTER

NAME: _____ DATE: _____
 ADDRESS: _____ AGE: _____
 ADDRESS: _____ WEARING GLASSES: Y N
 ADMINISTRATOR: _____ CONTACT LENSES: Y N

PRE-SCHOOL TEST

20/20, 20/30, 20/40, 20/50, 20/70, 20/100
 (Circle Acuity Level at which You are Testing) Must Read All Pictures Correctly
 TO WHICH PICTURE DO THE TABLE LEGS POINT?

RIGHT EYE TEST - FAR POINT

Left Occluder OFF

DIAL AT 02: BIRD
 DIAL AT 03: RABBIT
 DIAL AT 04: BOY
 DIAL AT 01: BINOCULARITY

PASS	FAIL

LEFT EYE TEST - FAR POINT

Right Occluder OFF

DIAL AT 04: BOY
 DIAL AT 03: GIRL
 DIAL AT 02: RABBIT

PASS	FAIL

	RIGHT EYE	20	20	20	20	20	20	20	20	20	20	20	20	20	20	20	20	20	20	LEFT EYE	20	20	20	20	20	20	20	20	20	20	20	20	20	20		
		200	100	70	50	40	30	20																												
PILOT LINE		E	m	E	Э	Э	W	M																												
20/50 BLUE	A	Э	E	W	Э	M	W																													
20/40 RED	B	W	W	Э	M	W	E																													
20/30 YELLOW	C	Э	M	E	Э	W																														
20/20 GREEN	D	E	W	E	Э	W	M																													

ELEMENTARY GRADE TESTS

ACUITY - FAR POINT

DIAL AT 05: RIGHT EYE
 DIAL AT 05: LEFT EYE

4 OUT OF 6 = PASS

Left Occluder OFF
 Right Occluder OFF

TEST FOR EXCESSIVE FARSIGHTEDNESS • PLUS LENS • FAR POINT (Insert Plus Lens Unit)

DIAL AT 05: LEFT EYE
 DIAL AT 05: RIGHT EYE

Right Occluder OFF
 Left Occluder OFF

(WITHDRAW PLUS LENS UNIT) A READING OF 20/30 OR 20/20 LINE THROUGH PLUS LENS = FAIL

MUSCLE BALANCE TEST - FAR POINT

DIAL AT 06: (1) Left Occluder OFF
 DIAL AT 06: (2) Both Occluders ON

(EXPLAIN SANDBOX)
 (RED BALL IN SANDBOX = PASS)

MUSCLE BALANCE TEST - NEAR POINT

DIAL AT 08: Both Occluders ON

(RED BALL STILL IN SANDBOX = PASS)

TEST RESULTS

20/
20/
PASS
FAIL
PASS
FAIL
PASS
FAIL

SECONDARY SCHOOL AND ADULT TEST

DIAL AT 09: ACUITY - FAR RIGHT EYE 20/____ LEFT EYE: 20/____ BOTH EYES: 20/____
 OPTIONAL A-1 2 B-5 C-2 6 D-6 E-1 6 F-BLANK

DIAL AT 07: COLOR PERCEPTION - FAR 1-R 2-L 3-B 4-T 5-B 6-L 7-T 8-R

DIAL AT 10: LATERAL PHORIA - FAR 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15

DIAL AT 11: VERTICAL PHORIA - FAR 1 2 3 4 5 6 7

DIAL AT 12: ACUITY - NEAR RIGHT EYE 20/____ LEFT EYE: 20/____ BOTH EYES: 20/____

DIAL AT 09: ACUITY - FAR	Score: 20/____		
LINE/ACUITY	LEFT EYE	BOTH EYES	RIGHT EYE
1 20/200	ZN	RO	HK
2 20/100	RKS	HNC	ZOD
3 20/70	HCDV	SKZO	RNDS
4 20/50	ZROD	NSCH	VZKN
5 20/40	KHSC	OZNR	DNVC
6 20/30	ONRZV	DKHCS	KDSON
7 20/20	SDCHN	VRZKO	HSNRD

DIAL AT 12: ACUITY - NEAR	Score: 20/____		
LINE/ACUITY	LEFT EYE	BOTH EYES	RIGHT EYE
1 20/100	SVC	NRK	HZO
2 20/70	RNZH	DOKV	CSZN
3 20/50	CKVD	SNZR	DOHC
4 20/40	VHRN	ODSK	NZCS
5 20/30	HSKRC	NZDOV	ZSHNK
6 20/20	ZONVR	HCKSD	VKCDSD
6 20/20	ZONVR	HCKSD	VKCDSD

Right Peripheral 85 ° 70 ° 55 ° NASAL 40 °
 Left Peripheral 85 ° 70 ° 55 ° NASAL 40 °



8600 W. Catalpa Avenue, Suite 703
 Chicago, IL 60656 - USA
 Phone: 1.800.344.9500
 Email: sales@stereo-optical.com
 Website : www.StereoOptical.com