# Optec<sup>®</sup> 2300 Armed Forces Vision Tester

# NSN 6540-01-375-9031



# **FEATURES INCLUDE:**

- Forehead Activator
- Forehead Activator Override
- Headrest Tissues
- Lens System For FAR and NEAR Testing
- Elevation Adjustment
- Observation Doors For Pointer Use & Test Slide Access
- FAR and NEAR Indicator Lights
- Test Dial To Change Slides In Viewing Area
- Carrying Handle

### The following tests can be performed using the vision tester.

# FAR TESTS

- Vertical Phoria
- Lateral Phoria
- Monocular Acuity (small letters)
- Monocular Acuity (large letters)
- Fusion and Depth Perception
- Alternate Monocular Acuity (small)
- Alternate Monocular Acuity (large)
- Alternate Fusion and Depth Perception
- Landolt Ring Acuity

# NEAR TESTS

- Vertical Phoria
- Lateral Phoria
- Monocular Acuity (small letters)
- Monocular Acuity (large letters)
- Binocular Acuity (small)
- Binocular Acuity (large)

**Accessories:** The unit comes with the following accessories: 50 record forms, 500 headrest tissues 1 dust cover, and 1 plastic laminated 81/2 by 11 scoring key. Includes two operators's manuals (which includes service manuals.)

**Warranty:** Free from defects in material and workmanship for a period of one year from date of purchase to original buyer. Guarantee does not cover light bulbs, mirrors and lenses, damage in transit, or damage caused by buyer's misuse, carelessness or neglect.

#### ALL U.S. MILITARY USE THE ABOVE ORDERING INFORMATION ALL NON U.S. MILITARY PLEASE CONTACT STEREO OPTICAL FOR ORDERING INFORMATION.



Stereo Optical Co., Inc.

773-867-0380 | 800-344-9500 | Fax: 773-867-0388 www.StereoOptical.com

Optec is a registered trademark of Stereo Optical Company. Optec 2300 manufactured on U.S. Patent #4,844,607 & #4,452,515