

STEREO OPTICAL REHAB G RECORD FORM FOR THE OPTEC[®] VISION TESTER

NAME: _____ DATE: _____
 ADDRESS: _____ AGE: _____
 ADDRESS: _____ WEARING GLASSES: Y N
 ADMINISTRATOR: _____ CONTACT LENSES: Y N

PERIPHERAL TEST L = 85° 70° 55° 45° NASAL 45° R = 85° 70° 55°


DIAL AT 01: F A C T FAR
Row A
Spatial
Frequency

DIAL AT 02: F A C T FAR
Row B
Spatial
Frequency

DIAL AT 03: F A C T FAR
Row C
Spatial
Frequency

DIAL AT 04: F A C T FAR
Row D
Spatial
Frequency

DIAL AT 05: F A C T FAR
Row E
Spatial
Frequency



DIAL AT 06: **ACUITY - FAR** Score: 20/ _____

LINE/ACUITY	LEFT EYE	BOTH EYES	RIGHT EYE
1 20/200	ZN	RO	HK
2 20/100	RKS	HNC	ZOD
3 20/70	HCDV	SKZO	RNDS
4 20/50	ZROD	NSCH	VZKN
5 20/40	KHSC	OZNR	DNVC
6 20/30	ONRZV	DKHCS	KDSON
7 20/20	SDCHN	VRZKO	HSNRD

DIAL AT 08: **ACUITY - NEAR** Score: 20/ _____

LINE/ACUITY	LEFT EYE	BOTH EYES	RIGHT EYE
1 20/100	SVC	NRK	HZO
2 20/70	RNZH	DOKV	CSZN
3 20/50	CKVD	SNZR	DOHC
4 20/40	VHRN	ODSK	NZCS
5 20/30	HSKRC	NZDOV	ZSHNK
6 20/20	ZONVR	HCSKD	VKCD S

DIAL AT 07: **COLOR PERCEPTION - FAR** A-1 2 B-5 C-2 6 D-6 E-1 6 F-BLANK PASS _____ FAIL _____

DIAL AT 09: **LATERAL PHORIA - FAR** 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15
 1 PRISM DIOPTERS ESOPHORIA ORTHO EXOPHORIA

DIAL AT 10: **DEPTH PERCEPTION - FAR** 1-B 2-L 3-B 4-T 5-T 6-L 7-R 8-L 9-R
 SECONDS OF ARC 400 200 100 70 50 40 30 25 20

- DIAL AT 11: **SIGN RECOGNITION / * DENOTES DEPTH PERCEPTION**
- | | | |
|---------------------------------------|-----------------------------|------------------------------------|
| 1. HILL | 5. REGULATION | 9. STOP SIGN |
| 2. HANDICAPPED PARKING | 6. * YIELD (323 Sec of Arc) | 10. RR CROSSING |
| 3. * DIVIDED HIGHWAY (581 Sec of Arc) | 7. NO U TURN | 11. * SCHOOL ZONE (145 Sec of Arc) |
| 4. NO PARKING | 8. DO NOT ENTER | 12. NO PASSING |

DIAL AT 12: **COLOR RECOGNITION**

LINE	1	2	3	4
LINE A	1=AMBER	2=GREEN	3=RED	4=BLUE
LINE B	1=RED	2=AMBER	3=BLUE	4=GREEN
LINE C	1=GREEN	2=BLUE	3=AMBER	4=RED
LINE D	1=BLUE	2=RED	3=GREEN	4=AMBER



1.800.344.9500 / 1.773.867.0380 / Fax: 1.773.867.0388
 sales@stereooptical.com / www.stereooptical.com
 8600 W. Catalpa Ave., Suite 703
 Chicago, IL 60656
 © Stereo Optical Co., Inc. 2017

REHAB G