



Optec® 1000P Vision Tester

Remote Control Vision Tester

Distance vision tests

Easily administered by ancillary office staff in less than 2 minutes

Features include:

- Test peripheral vision on the horizontal plane with stimuli at 85°, 70°, & 55° temporal & 45° nasal
- Offers a full range of tests including:
 - Acuity
 - Color
 - Muscle balance
- All test slides are produced on high-quality photographic film sealed between 2 glass plates to ensure quality images
- Requires less than 2 square feet of space and weighs less than 10 lbs (4.5 kg)
- Standard Optec 1000P package includes:
Adult and pediatric slides, plus lens (+1.75 D) and Allen test training flash card
- Additional optional accessories:
 - Carrying case
 - Plus lens (+2.25 D)
 - Intermediate lenses available in: 20 in, 22.5 in, 26 in, 31.5 in, & 39.5 in (50 cm, 57.15 cm, 66 cm, 80 cm, & 1 m)



-
- The Optec® 1000P with peripheral test and remote control includes: accessory kit*, and 1 standard slide package

*Accessory kit includes: dust cover, pointer, lens cleaner, headrest tissues, recording forms, and instruction manual
