

STEREO OPTICAL SPECIAL RECORD FORM FOR THE OPTEC® VISION TESTER

NAME: _____ DATE: _____
 ADDRESS: _____ AGE: _____
 ADDRESS: _____ WEARING GLASSES: Y N
 ADMINISTRATOR: _____ CONTACT LENSES: Y N

DIAL AT 01: ACUITY PATTI - FAR BINOCULAR

DIAL AT 02: ACUITY PATTI - FAR MONOCULAR

ACUITY	COLOR	BOTH EYES	ACUITY	COLOR	RIGHT EYE - LEFT EYE
20/100	WHITE	ST SQ H C A	20/100	WHITE	ST SQ C A H - H A C SQ ST
20/40	RED	SQ A ST H C	20/40	RED	C ST A H SQ - SQ H A ST C
20/32	AMBER	A H SQ C ST	20/32	AMBER	H A C SQ ST - ST SQ C A H
20/25	GREEN	SQ ST A H C	20/25	GREEN	SQ H A ST C - C ST A H SQ
20/20	BLUE	ST SQ C A H	20/20	BLUE	A SQ H C ST - ST C H SQ A

DIAL AT 03: ACUITY PATTI - NEAR BINOCULAR

DIAL AT 04: ACUITY PATTI - NEAR MONOCULAR

ACUITY	COLOR	BOTH EYES	ACUITY	COLOR	RIGHT EYE - LEFT EYE
20/100	WHITE	A ST C H SQ	20/100	WHITE	H ST A SQ C - C SQ A ST H
20/40	RED	ST SQ A C H	20/40	RED	A H SQ C ST - ST C SQ H A
20/32	AMBER	SQ C ST H A	20/32	AMBER	C SQ H ST A - A ST H SQ C
20/25	GREEN	ST A SQ C H	20/25	GREEN	ST C SQ A H - H A SQ C ST
20/20	BLUE	A ST H SQ C	20/20	BLUE	SQ ST C H A - A H C ST SQ

ST=Star, SQ=Square, C=Circle, A=Apple, H=House

DIAL AT 05: LEFT EYE Score: 20/_____ Right Occluder OFF

RIGHT EYE Score: 20/_____ Left Occluder OFF

TUMBLING "E" FAR	PILOT	20	20	20	20	20	20	PILOT	20	20	20	20	20	20	20	
		200	100	70	50	40	30		20	200	100	70	50	40	30	20
LINE	LINE	1	2	3	4	5	6	LINE	1	2	3	4	5	6	6	
LINE/ACUITY	COLOR	L	T	B	L	R	L	T	LINE/ACUITY	R	B	R	L	L	T	B
A 20/50	BLUE	R	T	R	L	T	B	A 20/50	L	R	T	L	B	T		
B 20/40	RED	B	R	B	L	L	R	B 20/40	T	T	L	B	T	R		
C 20/30	AMBER	L	T	T	L	R	B	C 20/30	L	B	R	R	L	T		
D 20/20	GREEN	R	L	T	L	B	T	D 20/20	R	T	R	L	T	B		

DIAL AT 06: MUSCLE BALANCE - FAR

Ball **IN** Box = **PASS**

Ball **OUT** of Box = **FAIL**

DIAL AT 07: COLOR PERCEPTION - FAR

1-R 2-L 3-B 4-T 5-B 6-L 7-T 8-R

DIAL AT 08: MUSCLE BALANCE - NEAR

Ball **IN** Box = **PASS**

Ball **OUT** of Box = **FAIL**

DIAL AT 09: ACUITY - FAR

Score: 20/_____

LINE/ACUITY	LEFT EYE	BOTH EYES	RIGHT EYE
1 20/200	Z N	R O	H K
2 20/100	R K S	H N C	Z O D
3 20/70	H C D V	S K Z O	R N D S
4 20/50	Z R O D	N S C H	V Z K N
5 20/40	K H S C	O Z N R	D N V C
6 20/30	O N R Z V	D K H C S	K D S O N
7 20/20	S D C H N	V R Z K O	H S N R D

DIAL AT 12: ACUITY - NEAR

Score: 20/_____

LINE/ACUITY	LEFT EYE	BOTH EYES	RIGHT EYE
1 20/100	S V C	N R K	H Z O
2 20/70	R N Z H	D O K V	C S Z N
3 20/50	C K V D	S N Z R	D O H C
4 20/40	V H R N	O D S K	N Z C S
5 20/30	H S K R C	N Z D O V	Z S H N K
6 20/20	Z O N V R	H C S K D	V K C D S

DIAL AT 10: LATERAL PHORIA - FAR

1 PRISM DIOPTERS

1	2	3	4	5	6	7	8	9	10	11	12	13	14	15
7	6	5	4	3	2	1	1	2	3	4	5	6	7	
ESOPHORIA						ORTHO			EXOPHORIA					

DIAL AT 11: VERTICAL PHORIA - FAR

1/2 PRISM DIOPTERS

1	2	3	4	5	6	7
1-1/2	1	1/2	0	1/2	1	1-1/2
LEFT HYPERPHORIA			ORTHO		RIGHT HYPERPHORIA	



**STEREO
OPTICAL**

1.800.344.9500 / 1.773.867.0380 / Fax: 1.773.867.0388
 sales@stereooptical.com / www.stereooptical.com
 8600 W Catalpa Avenue, Suite 703
 Chicago, IL 60656
 © Stereo Optical Co., Inc. 2017