

# Stereo Optical Slide Library



## 2000-001

### Far Acuity Binocular Landolt Rings


(20/200, 20/100, 20/70, 20/50, 20/40, 20/35, 20/30, 20/25, 20/22, 20/20, 20/18, 20/17, 20/15, 20/13)



## 2000-002

### Near Acuity Binocular Landolt Rings


(20/200, 20/100, 20/70, 20/50, 20/40, 20/35, 20/30, 20/25, 20/22, 20/20, 20/18, 20/17, 20/15, 20/13)



## 2000-003

### Far Acuity Binocular Letter

(20/200, 20/100, 20/70, 20/50, 20/40, 20/30, 20/20)



## 2000-004

### Far Acuity Tumbling E Monocular

(20/200, 20/100, 20/70, 20/50, 20/40, 20/30, 20/20)



## 2000-005

### Far Acuity Monocular Left Eye Landolt Rings

(20/200, 20/100, 20/70, 20/50, 20/40, 20/35, 20/30, 20/25, 20/22, 20/20, 20/18, 20/17, 20/15, 20/13)



## 2000-006

### Near Acuity Monocular Left Eye Landolt Rings

(20/200, 20/100, 20/70, 20/50, 20/40, 20/35, 20/30, 20/25, 20/22, 20/20, 20/18, 20/17, 20/15, 20/13)



**2000-007**  
Near Acuity Binocular Letter  
(20/100, 20/70, 20/50, 20/40, 20/30, 20/20)



**2000-008**  
Far Acuity Monocular Right Eye Landolt Rings  
(20/200, 20/100, 20/70, 20/50, 20/40, 20/35, 20/30, 20/25, 20/22, 20/20, 20/18, 20/17, 20/15, 20/13)



**2000-009**  
Near Acuity Monocular Right Eye Landolt Rings  
(20/200, 20/100, 20/70, 20/50, 20/40, 20/35, 20/30, 20/25, 20/22, 20/20, 20/18, 20/17, 20/15, 20/13)




**2000-010**  
Far Color Perception Binocular  
(Pseudo Ishihara)



**2000-011**  
Far Demonstration Mountains Binocular  
(Fixation for peripheral test)



**2000-012**  
Far Lateral Phoria Binocular  
(1 Diopter Increments)




**2000-013**  
Near Lateral Phoria Binocular  
(1-1/2 Diopter Increments)



**2000-014**  
Far Acuity Michigan Preschool Binocular  
20/70




**2000-015**  
Far Acuity Michigan Preschool Monocular  
20/30




**2000-016**  
Far Acuity Michigan Preschool Monocular  
20/30



**2000-017**  
Far Acuity Michigan Preschool Monocular  
20/30




**2000-018**  
Far Acuity Michigan Preschool Monocular  
20/30



**2000-019**  
Far Pediatric Color Perception Binocular  
Tumbling "E"




**2000-020**  
Far Muscle Balance Binocular  
(Combination Lateral and Vertical Phoria)



**2000-021**  
Far Multiple Phoria Binocular  
(2.3 Diopter Horizontal | 3.6 Diopter Vertical)



**2000-022**  
Near Muscle Balance Binocular  
(Combination Lateral and Vertical Phoria)



**2000-023**  
Near Multiple Phoria Binocular  
(2.7 Diopter Horizontal | 3.5 Diopter Vertical)



**2000-024**  
Far Stereo Depth Perception Binocular  
(400, 200, 100, 70, 50, 40, 30, 25, 20 Seconds of Arc)



**2000-025**  
Far Vertical Phoria Binocular  
(1/2 Diopter Increments)



**2000-026**  
Far Acuity HOTV Binocular  
20/70



**2000-027**  
Far Acuity HOTV Monocular  
20/40



**2000-028**  
Far Acuity HOTV Monocular  
20/40



**2000-029**  
Far Acuity HOTV Monocular  
20/40



**2000-030**  
Far Acuity HOTV Monocular  
20/30



**2000-031**  
Far Acuity HOTV Monocular  
20/30



**2000-032**  
Far Acuity HOTV Monocular  
20/30




**2000-033**  
Far Acuity HOTV Monocular  
20/30



**2000-034**  
Far Acuity HOTV Monocular  
20/30



**2000-035**  
Far Acuity Binocular Landolt Rings  
20/200, 20/15




**2000-042**  
Far Acuity Binocular  
(20/160, 20/125, 20/100, 20/80, 20/63, 20/50, 20/40, 20/32, 20/25, 20/20, 20/16, 20/12.5)




**2000-043**  
Far Astigmatic Clock Binocular




**2000-044**  
Far Acuity Binocular Number  
(20/200, 20/100, 20/70, 20/60, 20/50, 20/40, 20/30, 20/25, 20/20)




**2000-045**  
Near Acuity Binocular Number  
(20/200, 20/100, 20/70, 20/60, 20/50, 20/40, 20/30, 20/25, 20/20)



**2000-046**  
Far Acuity Binocular & Monocular Allen Test  
Binocular 20/100  
Monocular (20/50, 20/40, 20/30)



**2000-047**  
Far Acuity Binocular Letters  
(20/100, 20/70, 20/50, 20/40, 20/30)




**2000-048**  
Far Acuity Binocular Tumbling “E”  
(20/50, 20/40, 20/30, 20/20)




**2000-049**  
Far Acuity Binocular Letters  
(20/100, 20/70, 20/50, 20/40, 20/30, 20/20)




**2000-050**  
Far Acuity Binocular Letters  
(20/100, 20/75, 20/60, 20/50, 20/40, 20/20)




**2000-051**  
Far Acuity Binocular Letters  
(20/200, 20/100, 20/70, 20/50, 20/40, 20/30, 20/20)




**2000-052**  
Far Acuity Binocular Letters  
(20/100, 20/70, 20/50, 20/40, 20/30, 20/25, 20/20, 20/15)




**2000-054**  
Far Acuity Binocular Letters  
(20/100, 20/70, 20/50, 20/40, 20/30, 20/20)



**2000-055**  
Far Acuity Binocular Letter  
20/40 all lines




**2000-056**  
Far Acuity Binocular Letters  
Lines 1- 4 20/40  
Lines 5-6 20/100




**2000-057**  
Far Acuity Binocular Letters  
(20/200, 20/100, 20/70, 20/60, 20/50, 20/40, 20/30)



**2000-058**  
Far Acuity Binocular Tumbling "E"  
(20/100, 20/70, 20/60, 20/50, 20/40)




**2000-059**  
Far Acuity Binocular Letters  
(20/100, 20/70, 20/60, 20/50, 20/40, 20/30, 20/20)



**2000-060**  
Far Acuity Binocular Landolt C  
(20/100, 20/70, 20/60, 20/50, 20/40, 20/30, 20/20)




**2000-062**  
Far Acuity Monocular Numbers  
Top Row (20/200, 20/100, 20/70)  
(20/50, 20/40, 20/30, 20/25, 20/20, 20/15)



**2000-064**  
Far Acuity Binocular Numbers  
(20/70, 20/40, 20/30, 20/20)







**2000-065**  
**Far Acuity Binocular Numbers**  
 (20/200, 20/100, 20/70, 20/50, 20/40, 20/30, 20/20)



**2000-066**  
**Far Acuity Binocular Letters**  
 Top Row (20/200, 20/100)  
 (20/70, 20/60, 20/50, 20/40, 20/30, 20/25, 20/20)



**2000-067**  
**Far Acuity Binocular Numbers**  
 Top Row (20/200, 20/100)  
 (20/70, 20/60, 20/50, 20/40, 20/30, 20/25, 20/20)



**2000-068**  
**Far Acuity Monocular Numbers**  
 Top Row (20/200, 20/100, 20/70)  
 (20/50, 20/40, 20/30, 20/25, 20/20, 20/15)



**2000-070**  
**Far Acuity Monocular Tumbling "E"**  
 Top Row (20/200, 20/100)  
 (20/70, 20/50, 20/40, 20/30, 20/20)



**2000-071**  
**Near Acuity Binocular Tumbling "E"**  
 (20/50, 20/40, 20/30, 20/20)



**2000-072**  
**Near Acuity Allen Test**  
 Binocular 20/100  
 Monocular (20/50, 20/40, 20/30)



**2000-073**  
Near Acuity Binocular Letters  
(20/100, 20/70, 20/50, 20/40, 20/30, 20/20)



**2000-074**  
Near Acuity Monocular Tumbling "E"  
Top Row 20/200, 20/100, 20/70, 20/50, 20/40, 20/30, 20/20  
(20/50, 20/40, 20/30, 20/20)



**2000-075**  
Near Point Reading  
(.37,.50,.62,.75 M)



**2000-077**  
Far Point Color Perception (Pseudo Ishihara)  
Binocular Acuity Letter 20/40  
Binocular Acuity Landolt "C" 20/40  
Stereo Depth Perception (24 and 86 Seconds of Arc)



**2000-077C**  
Far Point color Perception (Pseudo Ishihara)  
Stereo Depth Perception (24 and 86 Seconds of Arc)




**2000-079**  
Far Acuity Binocular Letters 20/40  
Road Sign Recognition  
Color Recognition



**2000-080**  
Far Acuity Binocular Letters 20/40 Far Road  
Sign Recognition  
Far Color Perception (Pseudo Ishihara)



**2000-081**  
Far Acuity Binocular Letters 20/40  
Far Road Sign Recognition  
Far Color Perception (Pseudo Ishihara)



**2000-082**  
Far Color Perception Binocular  
(Pseudo Ishihara)




**2000-084**  
Far Color Recognition Binocular  
(Red, Green, Amber)



**2000-087**  
Far Fusion Binocular



**2000-088**  
Near Fusion Binocular



**2000-089**  
Far Acuity Binocular Landolt "C"  
(20/100, 20/75, 20/60, 20/50, 20/40, 20/20)



**2000-090**  
Far Acuity Binocular Michigan Preschool  
20/75



**2000-091**  
Near Acuity Binocular Michigan Preschool  
20/70



**2000-092**  
Far Acuity Monocular Michigan Preschool  
20/20




**2000-093**  
Far Acuity Monocular Michigan Preschool  
20/20



**2000-094**  
Far Acuity Binocular Michigan Preschool 20/20



**2000-095**  
Far Acuity Monocular Michigan Preschool  
20/20



**2000-096**  
Far Acuity Monocular Michigan Preschool  
20/40



**2000-097**  
Far Acuity Monocular Michigan Preschool  
20/40



**2000-098**

Far Acuity Monocular Michigan Preschool  
20/40




**2000-099**

Far Acuity Monocular Michigan Preschool  
20/40




**2000-100**

Far Acuity Monocular Michigan Preschool  
20/50



**2000-101**

Far Acuity Monocular Michigan Preschool  
20/50



**2000-102**

Far Acuity Monocular Michigan Preschool  
20/50



**2000-103**

Far Acuity Monocular Michigan Preschool  
20/50




**2000-104**

Far Acuity Monocular Michigan Preschool  
20/70



**2000-105**

Far Acuity Monocular Michigan Preschool  
20/70



**2000-106**

Far Acuity Binocular Michigan Preschool 20/70




**2000-107**

Far Acuity Monocular Michigan Preschool  
20/70



**2000-108**

Far Acuity Monocular Michigan Preschool  
20/100



**2000-109**

Far Acuity Monocular Michigan Preschool  
20/100




**2000-110**

Far Acuity Binocular Michigan Preschool  
20/100



**2000-111**

Far Acuity Monocular Michigan Preschool  
20/100



**2000-120**  
Near Stereo Depth Perception  
(400, 200, 100, 70, 50, 40, 30, 25, 20 Seconds of Arc)



**2000-125**  
Near Acuity Binocular Letters  
Top Row (20/200, 20/100, 20/70)  
(20/50, 20/40, 20/30, 20/25, 20/20, 20/15)



**2000-126**  
Far Acuity Monocular Letters  
Top Row (20/200, 20/100, 20/70)  
(20/50, 20/40, 20/30, 20/25, 20/20, 20/15)




**2000-128**  
Near Acuity Monocular Letters  
Top Row (20/200, 20/100, 20/70)  
(20/50, 20/40, 20/30, 20/25, 20/20, 20/15)



**2000-129**  
Far Acuity Monocular Letters  
Top Row (20/200, 20/100, 20/70)  
(20/50, 20/40, 20/30, 20/25, 20/20, 20/15)



**2000-132**  
Far Road Sign Recognition  
Depth Perception Binocular  
(620, 345 and 155 Seconds of Arc)



**2000-136**  
Far Acuity Low Vision Letters  
(20/206, 20/180, 20/163, 20/139, 20/127, 20/114)



**2000-137**  
 Far Acuity Low Vision "C"  
 (20/700, 20/600, 20/500, 20/400, 20/300, 20/200, 20/100)



**2000-138**  
 Far Acuity Low Vision Tumbling "E"  
 (20/800, 20/600, 20/400)



**2000-145**  
 Near Vertical Phoria Binocular  
 (1/2 Prism Diopter)



**2000-147**  
 Far Demonstration Highway Binocular  
 (Fixation for peripheral test)



**2000-148**  
 Far Road Sign Recognition and  
 Stereo Depth Perception Binocular  
 (620, 323 and 145 Seconds of Arc)



**2000-149**  
 Far Color Recognition  
 (Red, Yellow, Green, Blue)



**2000-182**  
 Far Visual Acuity Binocular HOTV  
 (20/100, 20/40, 20/32, 20/25, 20/20)





**2000-183**  
Far Visual Acuity Monocular HOTV  
(20/100, 20/40, 20/32, 20/25, 20/20)



**2000-184**  
Far Visual Acuity Binocular Numbers  
(20/40)



**2000-185**  
NEAR Acuity Monocular Letters  
(20/160, 20/125, 20/100, 20/80, 20/63, 20/50, 20/40, 20/32, 20/25, 20/20, 20/16, 20/12.5)




**2000-187**  
Near Visual Acuity Binocular HOTV  
(20/100, 20/40, 20/32, 20/25, 20/20)




**2000-188**  
Near Visual Acuity Monocular HOTV  
(20/100, 20/40, 20/32, 20/25, 20/20)



**2000-189**  
NEAR Acuity Binocular Letters  
(20/160, 20/125, 20/100, 20/80, 20/63, 20/50, 20/40, 20/32, 20/25, 20/20, 20/16, 20/12.5)



**2000-191**  
Far Binocular Contrast Letters  
25% ETDRS  
(20/160, 20/125, 20/100, 20/80, 20/63, 20/50, 20/40, 20/32, 20/25, 20/20, 20/16, 20/12.5)




## 2000-192

### Far Binocular Contrast Letters

#### 12.5% ETDRS

(20/160, 20/125, 20/100, 20/80, 20/63, 20/50 20/40, 20/32, 20/25, 20/20, 20/16, 20/12.5)




## 2000-193

### Far Binocular Contrast Letters

#### 50% ETDRS

(20/160, 20/125, 20/100, 20/80, 20/63, 20/50 20/40, 20/32, 20/25, 20/20, 20/16, 20/12.5)



## 2000-195

### NEAR Acuity Binocular Letters

#### 100% ETDRS


(20/160, 20/125, 20/100, 20/80, 20/64, 20/50 20/40, 20/32, 20/25, 20/20, 20/16, 20/12.5)



## 2000-253

### Far Vertical Phoria Binocular

#### (1/2 Diopter Increments)



## 2000-254

### Far Lateral Phoria Binocular

#### (1 Diopter Increments)



## 2000-255

### Far Acuity Binocular Checkerboard

(20/200, 20/100, 20/67, 20/50, 20/40, 20/33, 20/29, 20/25, 20/22, 20/20, 20/18, 20/17)



## 2000-256

### Far Acuity Monocular Right Eye Checkerboard

(20/200, 20/100, 20/67, 20/50, 20/40, 20/33, 20/29, 20/25, 20/22, 20/20, 20/18, 20/17)



**2000-257**  
 Far Acuity Monocular Left Eye Checkerboard  
 (20/200, 20/100, 20/67, 20/50, 20/40, 20/33, 20/29,  
 20/25, 20/22, 20/20, 20/18, 20/17)



**2000-258**  
 Far Stereo Depth Perception  
 (362, 83, 43, 32, 27, 19, 13, 12, 9.7 Seconds of Arc)




**2000-259**  
 Far Color Perception  
 (Pseudo Ishihara)



**2000-260**  
 Near Acuity Binocular Checkerboard  
 (20/200, 20/100, 20/67, 20/50, 20/40, 20/33, 20/29, 20/25, 20/22,  
 20/20, 20/18, 20/17)




**2000-261**  
 Near Acuity Monocular Right Eye  
 Checkerboard  
 (20/200, 20/100, 20/67, 20/50, 20/40, 20/33, 20/29,  
 20/25, 20/22, 20/20, 20/18, 20/17)



**2000-262**  
 Near Acuity Monocular Left Eye Checkerboard  
 (20/200, 20/100, 20/67, 20/50, 20/40, 20/33, 20/29,  
 20/25, 20/22, 20/20, 20/18, 20/17)



**2000-263**  
 Near Vertical Phoria Binocular  
 (1/2 Diopter Increments)




### 2000-264

Near Lateral Phoria Binocular  
(1 Diopter Increments)



### 2000-265

Far Acuity Monocular Tumbling "E"  
(20/200, 20/100, 20/50, 20/40, 20/30, 20/20)



### 2000-268

Far Acuity Binocular Letter  
(20/200, 20/100, 20/70, 20/50, 20/40, 20/30, 20/20)



### 3000-037

100% Far Acuity Monocular  
C.A.T.® Contrast Sensitivity  
(20/160, 20/125, 20/100, 20/80, 20/63, 20/50, 20/40, 20/32, 20/25, 20/20, 20/16, 20/12.5)



### 3000-038

50% Far Acuity Monocular  
C.A.T.® Contrast Sensitivity  
(20/160, 20/125, 20/100, 20/80, 20/63, 20/50, 20/40, 20/32, 20/25, 20/20, 20/16, 20/12.5)  
(Slides 037 to 041 sold as set only)



### 3000-039

25% Far Acuity Monocular  
C.A.T.® Contrast Sensitivity  
(20/160, 20/125, 20/100, 20/80, 20/63, 20/50, 20/40, 20/32, 20/25, 20/20, 20/16, 20/12.5)  
(Slides 037 to 041 sold as set only)



### 3000-040

12.5% Far Acuity Monocular  
C.A.T.® Contrast Sensitivity  
(20/160, 20/125, 20/100, 20/80, 20/63, 20/50, 20/40, 20/32, 20/25, 20/20, 20/16, 20/12.5)  
(Slides 037 to 041 sold as set only)




### 3000-041

6.25% Far Acuity Monocular

C.A.T.® Contrast Sensitivity

(20/160,20/125,20/100,20/80,20/63,20/50,20/40, 20/32,20/25,20/20,20/16, 20/12.5)

(Slides 037 to 041 sold as set only)




### 3000-171

Far F.A.C.T.® Contrast Sensitivity Monocular

1.5 Cycles Per Degree

(Slides 171 to 175 sold as set only)




### 3000-172

Far F.A.C.T. Contrast Sensitivity Monocular

3 Cycles Per Degree

(Slides 171 to 175 sold as set only)




### 3000-173

Far F.A.C.T. Contrast Sensitivity Monocular

6 Cycles Per Degree

(Slides 171 to 175 sold as set only)




### 3000-174

Far F.A.C.T. Contrast Sensitivity Monocular

12 Cycles Per Degree

(Slides 171 to 175 sold as set only)




### 3000-175

Far F.A.C.T.® Contrast Sensitivity Monocular

18 Cycles Per Degree

(Slides 171 to 175 sold as set only)




### 3000-176


Near F.A.C.T. Contrast Sensitivity Monocular

1.5 Cycles Per Degree

(Slides 176 to 178 sold as set only)



3000-177  
Near F.A.C.T. Contrast Sensitivity Monocular  
3 and 6 Cycles Per Degree  
*(Slides 176 to 178 sold as set only)*



3000-178  
Near F.A.C.T. Contrast Sensitivity Monocular  
12 & 18 Cycles Per Degree  
*(Slides 176 to 178 sold as set only)*



[info@stereooptical.com](mailto:info@stereooptical.com) | 773.867.0380 | [www.stereooptical.com](http://www.stereooptical.com)

© Stereo Optical Co. Inc. 2019. All rights reserved.