




Optec Rehab Slide Package




Slide #1 (3000-171)** Far Point
Functional Acuity Contrast Test
1.5 Cycles Per Degree Monocular




Slide #2 (3000-172) Far Point
Functional Acuity Contrast Test
3 Cycles Per Degree Monocular




Slide #3 (3000-173) Far Point
Functional Acuity Contrast Test
6 Cycles Per Degree Monocular




Slide #4 (3000-174) Far Point
Functional Acuity Contrast Test
12 Cycles Per Degree Monocular




Slide #5 (3000-175) Far Point
Functional Acuity Contrast Test
18 Cycles Per Degree Monocular




Slide #6 (2000-003) Far Point
Acuity Right Eye & Left Eye Monocular
or Both Eyes Simultaneously (20/200-20/20)




Slide #7 (2000-010) Far Point
Color Perception (Pseudo Ishihara)



Slide #8 (2000-007) Near Point
Acuity Right Eye & Left Eye Monocular
or Both Eyes Simultaneously (20/100-20/20)




Slide #9 (2000-012) Far Point
Lateral Phoria (1 Diopter Increments)



Slide #10 (2000-024) Far Point
stereo Depth Perception
(400-20 Seconds of Arc)



Slide #11 (2000-148) Far Point
Road Sign Recognition & Depth Perception
(620-323-145 Seconds of Arc)



Slide #12 (2000-149) Far Point
Color Recognition (Red, Yellow, Green, Blue)

** Enlarged to show pattern