



Stereo Optical Patti Pic Slide Package


For the complete vision screening needs from preschool through adult. Ideal for School, Athletic Physicals, Employment Physicals, and General Check-Ups.




Slide #1 (2000-182P)
Far Acuity Binocular Patti Pic
(20/100, 20/40, 20/32, 20/25, 20/20)




Slide #2 (2000-183P)
Far Acuity Monocular Patti Pic
(20/100, 20/40, 20/32, 20/25, 20/20)




Slide #3 (2000-187P)
Near Acuity Binocular Patti Pic
(20/100, 20/40, 20/32, 20/25, 20/20)




Slide #4 (2000-188P)
Near Acuity Monocular Patti Pic
(20/100, 20/40, 20/32, 20/25, 20/20)




Slide #5 (2000-004) Far Point
Acuity Right Eye & Left Eye Monocular,
Right Eye & Left Eye Plus Lens Test (20/200-20/20)




Slide #6 (2000-020) Far Point Muscle Balance
(Combination Lateral & Vertical Phoria)



Slide #7 (2000-019) Far Point
Color Perception (Tumbling "E")



Slide #8 (2000-022) Near Point
Muscle Balance
(Combination Lateral & Vertical Phoria)




Slide #9 (2000-003) Far Point
Acuity Right Eye & Left Eye Monocular
or Both Eyes Simultaneously (20/200-20/20)



Slide #10 (2000-012) Far Point
Lateral Phoria (1 Diopter Increments)



Slide #11 (2000-025) Far Point Vertical
Phoria (1/2 Diopter Increments)



Slide #12 (2000-007) Near Point
Acuity Right Eye & Left Eye Monocular
or Both Eyes Simultaneously (20/100-20/20)