

# STEREO OPTICAL DRIVER REHAB RECORD FORM FOR THE OPTEC® VISION TESTER

NAME: \_\_\_\_\_ DATE: \_\_\_\_\_  
 ADDRESS: \_\_\_\_\_ AGE: \_\_\_\_\_  
 CITY, STATE, ZIP: \_\_\_\_\_ WEARING GLASSES: Y N  
 Administrator: \_\_\_\_\_ CONTACT LENSES: Y N

**PERIPHERAL TEST** L = 85° 70° 55° 45°

**DIAL AT 01: ACUITY – FAR**

	Score: 20/ ____	Score: 20/ ____	Score: 20/ ____
LINE/ACUITY	LEFT EYE	BOTH EYES	RIGHT EYE
1 20/200	Z N	R O	H K
2 20/100	R K S	H N C	Z O D
3 20/70	H C D V	S K Z O	R N D S
4 20/50	Z R O D	N S C H	V Z K N
5 20/40	K H S C	O Z N R	D N V C
6 20/30	O N R Z V	D K H C S	K D S O N
7 20/20	S D C H N	V R Z K O	H S N R D

R = 85° 70° 55° 45°

**DIAL AT 03: ACUITY – NEAR**

	Score: 20/ ____	Score: 20/ ____	Score: 20/ ____
LINE/ACUITY	LEFT EYE	BOTH EYES	RIGHT EYE
1 20/100	S V C	N R K	H Z O
2 20/70	R N Z H	D O K V	C S Z N
3 20/50	C K V D	S N Z R	D O H C
4 20/40	V H R N	O D S K	N Z C S
5 20/30	H S K R C	N Z D O V	Z S H N K
6 20/20	Z O N V R	H C S K D	V K C D S

**DIAL AT 02: COLOR PERCEPTION – FAR** A-1 2 B-5 C-2 6 D-6 E-1 6 F-BLANK PASS FAIL

**DIAL AT 04: LATERAL PHORIA – FAR** Score: \_\_\_\_

1	2	3	4	5	6	7	8	9	10	11	12	13	14	15
ESOPHORIA						ORTHO						EXOPHORIA		

**DIAL AT 05: VERTICAL PHORIA – FAR** Score: \_\_\_\_

1	2	3	4	5	6	7
LEFT HYPERPHORIA			ORTHO		RIGHT HYPERPHORIA	

**DIAL AT 06: LATERAL PHORIA – NEAR** Score: \_\_\_\_

1	2	3	4	5	6	7	8	9	10	11	12	13	14	15
ESOPHORIA						ORTHO						EXOPHORIA		

**DIAL AT 07: DEPTH PERCEPTION – FAR** SECONDS OF ARC

1-B	2-L	3-B	4-T	5-T	6-L	7-R	8-L	9-R	Score:
400	200	100	70	50	40	30	25	20	____

**DIAL AT 08: FUSION – FAR** 3 cubes BOTH eyes = **PASS** Other than 3 cubes = **FAIL** PASS FAIL

**DIAL AT 09: TUMBLING “E” FAR**

LINE/ACUITY	COLOR	LEFT EYE Score: 20/ ____ Right Occluder OFF								RIGHT EYE Score: 20/ ____ Left Occluder OFF														
		PILOT		20	20	20	20	20	20	20	20	PILOT		20	20	20	20	20	20	20	20			
		200	100	70	50	40	30	20	200	100	70	50	40	30	20	200	100	70	50	40	30	20		
LINE/ACUITY	COLOR	L	T	B	L	R	L	T	LINE/ACUITY	R	B	R	L	L	T	B	LINE/ACUITY	R	B	R	L	L	T	B
A 20/50	BLUE								A 20/50							A 20/50								
B 20/40	RED								B 20/40							B 20/40								
C 20/30	AMBER								C 20/30							C 20/30								
D 20/20	GREEN								D 20/20							D 20/20								

**DIAL AT 10: ACUITY – FAR**

	20/ ____	20/ ____	20/ ____	
LINE/ACUITY	COLOR	LEFT	BOTH	RIGHT
1 20/200	RED	3 9	6 5	2 8
2 20/100	GREEN	8 5 2	3 9 5	6 2 8
3 20/70	RED	5 3 6 2	9 8 5 6	2 3 6 5
4 20/50	AMBER	9 2 8 6	5 9 6 3	2 8 5 9
5 20/40	GREEN	5 8 6 9	2 6 8 3	9 5 3 2
6 20/30	RED	3 2 5 9 6	6 8 3 5 2	8 2 6 3 5
7 20/20	AMBER	2 9 8 6 5	3 5 2 6 9	6 5 3 8 2

**DIAL AT 11: SIGN RECOGNITION AND DEPTH PERCEPTION**

<b>LINE 1</b>	Y	N
1. HILL		
2. HANDICAPPED PARKING		
3. * DIVIDED HIGHWAY (581 Sec of Arc)		
4. NO PARKING		
<b>LINE 2</b>	Y	N
5. REGULATION		
6. * YIELD (323 Sec of Arc)		
7. NO U TURN		
8. DO NOT ENTER		
<b>LINE 3</b>	Y	N
9. STOP SIGN		
10. RR CROSSING		
11. * SCHOOL ZONE (145 Sec of Arc)		
12. NO PASSING		

**DIAL AT 12: COLOR RECOGNITION**

LINE	1=AMBER	2=GREEN	3=RED	4=BLUE	Y	N
LINE A						
LINE B						
LINE C						
LINE D						



**STEREO OPTICAL CO., INC.**

8600 W. Catalpa Avenue, Suite 703  
 Chicago, IL 60656  
 773.867.0380

www.stereooptical.com | sales@stereooptical.com

© Stereo Optical Co., Inc. 10/2020. All rights reserved. Unless indicated otherwise, all trademarks are the property of Stereo Optical Co., Inc. in the United States and other countries.