



The Optec® 5000G INT Vision Screener

Perform comprehensive visual screening in minimal time. Screen under a variety of conditions, including monocular, binocular, near, distance, and intermediate. Test for adult and pediatric visual acuities, color and depth perception, as well as vertical and lateral phorias. The Optec 5000G INT vision screener can also function as a peripheral vision test.



FRONT VIEW

FEATURES INCLUDE:

- Quick, repeatable and accurate
- Manual control dial and use of color-corrected LED's
- Test with or without glare (non-calibrated)
- Lightweight and portable
- See and point observation doors on both sides
- No training or certification needed

*Includes a slide package with a set of 12 test slides as well as the option to create a custom slide package from a library of more than 150 test slides**



BACK VIEW

STANDARD TEST PACKAGES HAVE BEEN DESIGNED FOR:

- Ophthalmology
- Public Health
- Optometry
- Primary Care
- Occupational Medicine
- D.O.T. (Dept. of Transportation)
- Driver Rehabilitation
- Schools
- Pediatrics
- Sports/Athletics
- F.A.A. (Federal Aviation Administration)
- Customizable Test Packages

Ancillary staff can test Acuities (Near/Far/Monocular/Binocular), Ishihara Color, Stereopsis (Depth Perception), Phorias, Astigmatic Clock, and much more in under 3 minutes!

Lightweight and portable, weighing less than 15 lbs (6.8 kg), requiring only 2 square ft (0.6 m) of space, has scratch resistant and washable surface made of ABS plastic.

*2-year warranty

**The F.A.C.T. (contrast sensitivity) test slides are not available with Optec 5000G INT.



- | | |
|---------------------------|--|
| 1. Forehead Activator | 7. Far and Near Indicators |
| 2. Headrest Tissue Holder | 8. Test Dial and Knob |
| 3. Lens System | 9. Receptacle for Power Cord and Control Panel |
| 4. Instrument Base | 10. Instrument Body |
| 5. Elevation Adjustment | 11. Carrying Handle |
| 6. Observation Door | 12. Ancillary Lens Holder |

OPTEC® 5000G INT WITH GLARE INCLUDES:

- Accessory Kit*
- One standard slide package
- Peripheral test

*Accessory kit includes: dust cover, pointer, lens cleaner, headrest tissues, recording forms, and instruction manual