


STEREO OPTICAL RECORD FORM FOR THE OPTEC® VISION TESTER

NAME: _____ DATE: _____
 ADDRESS: _____ AGE: _____
 CITY, STATE, ZIP: _____ WEARING GLASSES: Y N
 ADMINISTRATOR: _____ CONTACT LENSES: Y N

PERIPHERAL TEST L = 85° 70° 55° 45° R = 85° 70° 55° 45°

- DIAL AT 01: F A C T F A R
Row A
Spatial
Frequency
- DIAL AT 02: F A C T F A R
Row B
Spatial
Frequency
- DIAL AT 03: F A C T F A R
Row C
Spatial
Frequency
- DIAL AT 04: F A C T F A R
Row D
Spatial
Frequency
- DIAL AT 05: F A C T F A R
Row E
Spatial
Frequency



Day Night (button).
Glare No Glare

Complete page 2 if using more than one test lighting condition.

DIAL AT 06: ACUITY – FAR

	Score: 20/___	Score: 20/___	Score: 20/___
LINE/ACUITY	LEFT EYE	BOTH EYES	RIGHT EYE
1 20/200	Z N	R O	H K
2 20/100	R K S	H N C	Z O D
3 20/70	H C D V	S K Z O	R N D S
4 20/50	Z R O D	N S C H	V Z K N
5 20/40	K H S C	O Z N R	D N V C
6 20/30	O N R Z V	D K H C S	K D S O N
7 20/20	S D C H N	V R Z K O	H S N R D

DIAL AT 08: ACUITY – NEAR

	Score: 20/___	Score: 20/___	Score: 20/___
LINE/ACUITY	LEFT EYE	BOTH EYES	RIGHT EYE
1 20/100	S V C	N R K	H Z O
2 20/70	R N Z H	D O K V	C S Z N
3 20/50	C K V D	S N Z R	D O H C
4 20/40	V H R N	O D S K	N Z C S
5 20/30	H S K R C	N Z D O V	Z S H N K
6 20/20	Z O N V R	H C S K D	V K C D S

DIAL AT 07: COLOR PERCEPTION – FAR A-1 2 B-5 C-2 6 D-6 E-1 6 F-BLANK PASS FAIL

DIAL AT 09: LATERAL PHORIA – FAR Score: ___

<input type="checkbox"/> 1	<input type="checkbox"/> 2	<input type="checkbox"/> 3	4	5	6	7	8	9	10	11	12	<input type="checkbox"/> 13	<input type="checkbox"/> 14	<input type="checkbox"/> 15
ESOPHORIA						ORTHO						EXOPHORIA		

PRISM DIOPTERS

DIAL AT 10: DEPTH PERCEPTION – FAR

SECONDS OF ARC	<u>1-B</u>	<u>2-L</u>	<u>3-B</u>	<u>4-T</u>	<u>5-T</u>	<u>6-L</u>	<u>7-R</u>	<u>8-L</u>	<u>9-R</u>	Score:
	400	200	100	70	50	40	30	25	20	___

DIAL AT 11: SIGN RECOGNITION / * DENOTES DEPTH PERCEPTION

LINE 1	Y	N	LINE 2	Y	N	LINE 3	Y	N
1. HILL			5. REGULATION			9. STOP SIGN		
2. HANDICAPPED PARKING			6. *YIELD (323 Sec of Arc)			10. RR CROSSING		
3. * DIVIDED HIGHWAY (581 Sec of Arc)			7. NO U TURN			11. * SCHOOL ZONE (145 Sec of Arc)		
4. NO PARKING			8. DO NOT ENTER			12. NO PASSING		

DIAL AT 12: COLOR RECOGNITION

LINE A	1=AMBER	2=GREEN	3=RED	4=BLUE	Y	N
LINE B	1=RED	2=AMBER	3=BLUE	4=GREEN	Y	N
LINE C	1=GREEN	2=BLUE	3=AMBER	4=RED	Y	N
LINE D	1=BLUE	2=RED	3=GREEN	4=AMBER	Y	N



STEREO OPTICAL


STEREO OPTICAL CO., INC.

8600 W. Catalpa Avenue, Suite 703
Chicago, IL 60656
773.867.0380


www.stereooptical.com | info@stereooptical.com

© Stereo Optical Co., Inc. 10/2020. All rights reserved. Unless indicated otherwise, all trademarks are the property of Stereo Optical Co., Inc. in the United States and other countries.


NIGHT TESTING WITHOUT GLARE




NIGHT TESTING WITH GLARE



DAY TESTING WITHOUT GLARE



DAY TESTING WITH GLARE



STEREO OPTICAL CO., INC.

8600 W. Catalpa Avenue, Suite 703

Chicago, IL 60656

773.867.0380

www.stereooptical.com | info@stereooptical.com

© Stereo Optical Co., Inc. 09/2020. All rights reserved. Unless indicated otherwise, all trademarks are the property of Stereo Optical Co., Inc. in the United States and other countries.