

STEREO OPTICAL DRIVER REHAB RECORD FORM FOR THE OPTEC® VISION TESTER

NAME: _____ DATE: _____
 ADDRESS: _____ AGE: _____
 ADDRESS: _____ WEARING GLASSES: Y N
 ADMINISTRATOR: _____ CONTACT LENSES: Y N

PERIPHERAL TEST $L = 85^\circ 70^\circ 55^\circ$ 45° NASAL 45° $R = 85^\circ 70^\circ 55^\circ$

DIAL AT 01: ACUITY - FAR Score: 20/_____

LINE/ACUITY	LEFT EYE	BOTH EYES	RIGHT EYE
1 20/200	Z N	R O	H K
2 20/100	R K S	H N C	Z O D
3 20/70	H C D V	S K Z O	R N D S
4 20/50	Z R O D	N S C H	V Z K N
5 20/40	K H S C	O Z N R	D N V C
6 20/30	O N R Z V	D K H C S	K D S O N
7 20/20	S D C H N	V R Z K O	H S N R D

DIAL AT 03: ACUITY - NEAR Score: 20/_____

LINE/ACUITY	LEFT EYE	BOTH EYES	RIGHT EYE
1 20/100	S V C	N R K	H Z O
2 20/70	R N Z H	D O K V	C S Z N
3 20/50	C K V D	S N Z R	D O H C
4 20/40	V H R N	O D S K	N Z C S
5 20/30	H S K R C	N Z D O V	Z S H N K
6 20/20	Z O N V R	H C S K D	V K C D S

DIAL AT 02: COLOR PERCEPTION - FAR A-1 2 B-5 C-2 6 D-6 E-1 6 F-BLANK PASS ____ FAIL ____

DIAL AT 04: LATERAL PHORIA - FAR 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15
 1 PRISM DIOPTERS E S O P H O R I A O R T H O E X O P H O R I A

DIAL AT 05: VERTICAL PHORIA - FAR 1 2 3 4 5 6 7
 1/2 PRISM DIOPTERS L E F T H Y P E R P H O R I A O R T H O R I G H T H Y P E R P H O R I A

DIAL AT 06: LATERAL PHORIA - NEAR 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15
 1 PRISM DIOPTERS E S O P H O R I A O R T H O E X O P H O R I A

DIAL AT 07: DEPTH PERCEPTION - FAR 1-B 2-L 3-B 4-T 5-T 6-L 7-R 8-L 9-R
 SECONDS OF ARC 400 200 100 70 50 40 30 25 20

DIAL AT 08: FUSION - FAR 3 cubes BOTH eyes = PASS Other than 3 cubes = FAIL

DIAL AT 09: TUMBLING "E" FAR LEFT EYE Score: 20/_____ Right Occluder OFF RIGHT EYE Score: 20/_____ Left Occluder OFF

LINE/ACUITY	COLOR	LEFT EYE								RIGHT EYE							
		L	T	B	L	R	L	T	B	R	B	R	L	L	T	B	
A 20/50	BLUE	200	100	70	50	40	30	20	200	100	70	50	40	30	20		
B 20/40	RED																
C 20/30	AMBER																
D 20/20	GREEN																

DIAL AT 10: ACUITY - FAR Score: 20/_____

LINE/ACUITY	COLOR	LEFT	BOTH	RIGHT
1 20/200	RED	39	65	28
2 20/100	GREEN	8 5 2	3 9 5	6 2 8
3 20/70	RED	5 3 6 2	9 8 5 6	2 3 6 5
4 20/50	AMBER	9 2 8 6	5 9 6 3	2 8 5 9
5 20/40	GREEN	5 8 6 9	2 6 8 3	9 5 3 2
6 20/30	RED	3 2 5 9 6	6 8 3 5 2	8 2 6 3 5
7 20/20	AMBER	2 9 8 6 5	3 5 2 6 9	6 5 3 8 2

- DIAL AT 11: SIGN RECOGNITION AND DEPTH PERCEPTION**
1. HILL
 2. HANDICAPPED PARKING
 3. * DIVIDED HIGHWAY (581 Sec of Arc)
 4. NO PARKING
 5. REGULATION
 6. * YIELD (323 Sec of Arc)
 7. NO U TURN
 8. DO NOT ENTER
 9. STOP SIGN
 10. RR CROSSING
 11. * SCHOOL ZONE (145 Sec of Arc)
 12. NO PASSING

DIAL AT 12: COLOR RECOGNITION

LINE A	1=AMBER	2=GREEN	3=RED	4=BLUE
LINE B	1=RED	2=AMBER	3=BLUE	4=GREEN
LINE C	1=GREEN	2=BLUE	3=AMBER	4=RED
LINE D	1=BLUE	2=RED	3=GREEN	4=AMBER

© Stereo Optical Co., Inc. 05/2019. All rights reserved. Unless indicated otherwise, all trademarks are the property of Stereo Optical Co., Inc. in the United States and other countries.